1. MANUEL D'UTILISATION 2. USER GUIDE 3. GEBRAUCHSANLEITUNG 4. MANUAL DE USO 5. ISTRUZIONI D'USO 6. INSTRUÇÕES DE USO  
7. BRUKSANVISNING 8. GEBRUIKSAANWIJZING 9. INSTRUKCJA OBSŁUGI 10. ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ 11. ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ  
12. UŽIVATELSKÝ MANUÁL 13. 用戶手冊 14. 使用说明 15. 使用説明書


# Respireo **SOFT** Facial

1. MASQUE BUCCO-NASAL 2. FACIAL MASK 3. MUND-NASEN-MASKE 4. MASCARILLA ORO-NASAL 5. MASCHERA ORO-NASALE  
6. MÁSCARA ORO-NASAL 7. HELMmask 8. MOND-NEUSMASKER 9. MASKA USTNO-NOSOWA 10. ΠΙΝΙΚΗ ΣΤΟΜΑΤΙΚΗ ΜΑΣΚΑ  
11. ПОТО-НОСОВАЯ МАСКА 12. NOSO-ÚSTNÍ MASKA 13. 口鼻面膜 14. 口鼻面罩 15. 鼻口腔マスク




1. FRANÇAIS.....	p. 17
2. ENGLISH.....	p. 20
3. DEUTSCH.....	p. 23
4. ESPAÑOL.....	p. 27
5. ITALIANO .....	p. 30
6. PORTUGUÊS.....	p. 33
7. SVENSKA .....	p. 36
8. NEDERLANDS.....	p. 39
9. POLSKI.....	p. 42
10. ΕΛΛΗΝΙΚΑ.....	p. 46
11. РУССКИЙ.....	p. 49
12. ČEŠTINA.....	p. 53
13. 中國傳統.....	p. 56
14. 中文 .....	p. 59
15. 日本語.....	p. 62

**A**


**B**

1. COMMENT PORTER LE MASQUE 2. HOW TO WEAR THE MASK 3. KORREKTES ANLEGEN DER MASKE 4. CÓMO UTILIZAR LA MÁSCARA  
5. COME INDOSSARE LA MASCHERA 6. COME VESTIR A MÁSCARA 7. TA PÅ MASKEN SÅ HÅR 8. HET MASKER OPZETTEN 9. JAK ZAŁOŻYĆ  
MASKĘ 10. ΠΩΣ ΦΟΡΑΜΕ ΤΗ ΜΑΣΚΑ 11. КАК НАДЕВАТЬ МАСКУ 12. JAK MASKU NASADIT 13. 如何佩戴鼻罩 14. 如何佩戴鼻罩 15. マスクの着用方法

**1****2****3**

**B**


4

5


6





1. COMMENT AJUSTER LE MASQUE 2. HOW TO ADJUST THE MASK 3. KORREKTES EINSTELLEN DER MASKE 4. CÓMO REGULAR LA MÁSCARA 5. COME REGOLARE LA MASCHERA 6. COME REGULAR A MÁSCARA 7. JUSTERA MASKEN SÅ HÄR 8. HET MASKER REGELEN 9. JAK WYREGULOWAĆ MASKĘ 10. ΠΩΣ ΡΥΘΜΙΖΟΥΜΕ ΤΗ ΜΑΣΚΑ 11. ΚΑΚ ΡΕΓΥΛΙΡΟΒΑΤΉ ΜΑΣΚΥ 12. JAK MASKU UPRAVIT REGELEN 13. 如何調整鼻罩 14. 如何調整鼻罩 15. マスクの調整方法


1. Remarque : ne serrez pas trop les branches 2. Note: do not over-tighten the straps 3. Hinweis: Ziehen Sie die Gurte nicht zu fest an 4. Nota: no ajustar excesivamente las correas 5. Nota: non stringere eccessivamente le cinghie 6. Nota: não apertar excessivamente as correias 7. Anmärkning: dra inte åt remmarna för hårt 8. Opmerking: de riempjes niet te veel aanspannen 9. Uwaga: nie należy nadmierne dociskać pasków 10. Σημείωση: μην σφίγγετε υπερβολικά τις ζώνες 11. Примечание: не затягивайте чрезмерно ремни 12. Poznámka: popruhy příliš neutahujte 13. 備註: 勿過度收緊繫帶 14. 注: 不要將皮帶拉得過緊 15. 備考: ベルトの締め付け過ぎに注意してください。

**D**



1. COMMENT SE PRÉPARER À LA THÉRAPIE 2. HOW TO PREPARE FOR TREATMENT 3. VORBEREITUNG AUF DIE THERAPIE 4. CÓMO PREPARARSE PARA LA TERAPIA 5. COME PREPARARSI ALLA TERAPIA 6. COMO PREPARAR-SE À TERAPIA 7. FÖRBERED DIG FÖR TERAPIN SÅ HÄR 8. ZICH VOORBEREIDEN OP DE THERAPIE 9. JAK PRZYKOTOWAĆ SIĘ DO TERAPII 10. ΠΡΟΕΤΟΙΜΑΣΙΕΣ ΘΕΡΑΠΕΙΑΣ 11. КАК ПОДГОТОВИТЬСЯ К ТЕРАПИИ 12. JAK SE PŘIPRAVIT NA TERAPII 13. 如何準備接受治療 14. 如何准备治疗 15. 治療への準備方法

4

5


**1. Remarques 2. Notes 3. Hinweise 4. Notas 5. Note 6. Notas 7. Anmärkning 8. Opmerkingen 9. Uwagi 10. Σημειώσεις  
11. Примечания 12. Poznámky 13. 注意 14. 注意 15. 注意事項**

- ▶ **1.** En cas de plis sur la bulle (a), répétez l'opération **2.** If the cushion (a) is wrinkled, repeat the operation **3.** Wenn das Kissen (a) Falten aufweist, wiederholen Sie den Vorgang **4.** Si el cojín (a) presenta arrugas, repita la operación **5.** Se il cuscino (a) presenta grinze, ripetere l'operazione **6.** Se a almofada (a) estiver enrugada, repetir a operação **7.** Om dyndan (a) uppvisar skrynklor, ska du upprepa ingreppet **8.** Als plooiën op het kussen (a) aanwezig zijn, de handling herhalen **9.** Jeśli poduszka (a) marszczy się, powtórzyć czynność **10.** Αν το μαξιλαράκι (a) παρουσιάζει ζάρες, επαναλάβετε τη διαδικασία **11.** Если на подушке (a) имеются складки (a), повторите операцию **12.** Pokud polštářek (a) vykazuje záhyby, operaci opakujte **13.** 如果墊子 (a) 上有褶皺，請重複操作 **14.** 如衬垫(a) 有皱纹，請重複操作 **15.** クッション(a) にしわが寄っている場合は、この作業を繰り返してください
- ▶ **1.** En cas de fuite, resserrez légèrement les branches du harnais (e) **2.** If there are leaks, tighten the headgear straps (e) **3.** Wenn es zu Undichtigkeiten kommt, ziehen Sie die Gurte des Kopfteils leicht an(e) **4.** En caso de pérdidas, ajuste ligeramente las correas de la cofia (e) **5.** In caso di perdite, stringere leggermente le cinghie della cuffia (e) **6.** Em caso de vazamentos, apertar ligeiramente as correias do protetor de cabeça (e) **7.** Vid läckage, dra åt hättans remmar något (e) **8.** Bij verlies de riempjes van de kap (e) enigszins aanspannen **9.** W przypadku wycieków, lekko docisnąć pasy uprząży (e) **10.** Σε περίπτωση διαρροών, σφίξτε ελαφρώς τις ζώνες του ακουστικού (e) **11.** в случае утечек слегка затяните ремни пращи (e) **12.** V případě úniků jemně přitlačte popruhy krytu (e) **13.** 如果出現洩漏，輕輕地收緊頭套 (e) 的繫帶 **14.** 如有漏气现象，请轻轻拉紧鼻罩的皮带(e) **15.** 漏れが見つかった場合は、ヘッドギア(e)のベルトを軽く締め付けてください


1. COMMENT ENLEVER ET DÉMONTÉ LE MASQUE 2. HOW TO REMOVE AND DISASSEMBLE THE MASK 3. KORREKTES ABNEHMEN UND DEMONTIEREN DER MASKE 4. CÓMO EXTRAER Y DESMONTAR LA MÁSCARA 5. COME RIMUOVERE E SMONTARE LA MASCHERA 6. COMO REMOVER E DESMONTAR A MÁSCARA 7. TA AV OCH MONTERA NED MASKEN SÅ HÄR 8. HET MASKER VERWIJDEREN EN DEMONTEREN 9. JAK ZDJAĆ I ZDEMONTOWAĆ MASKĘ 10. ΠΩΣ ΑΦΑΙΡΟΥΜΕ ΚΑΙ ΑΠΟΣΥΝΑΡΜΟΛΟΓΟΥΜΕ ΤΗ ΜΑΣΚΑ 11. ΚΑΚ ΣΗΜΑΤΩ Ι ΡΑΖΒΙΡΑΤΨ ΜΑΣΚΥ 12. JAK MASKU ODSTRANIT A DEMONTOVA 13. 如何取下和拆卸鼻罩 14. 如何移除和拆卸鼻罩 15. マスクの取り外しと解体方法


1. Remarque: ne retirez pas le masque sans avoir relâché les branches inférieures du harnais (e). 2. Note: do not remove the mask without first disconnecting the lower headgear straps (e). 3. Hinweis: nehmen Sie die Maske nicht ab, ohne die unteren Kopfteil-Gurte zu lösen. (e). 4. Nota: no extraiga la máscara sin haber desconectado las correas inferiores de la cofia (e). 5. Nota: non rimuovere la maschera senza aver scollato le cinghie inferiori della cuffia (e). 6. Nota: não remover a máscara sem ter desencalhado as correias inferiores do protetor de cabeça (e). 7. Anmärkning: ta inte av masken innan du har avlägsnat de nedre remmarna från hättan (e). 8. Opmerking: het masker niet verwijderen zonder de onderste riempjes van de kap (e) los te maken. 9. Uwaga: nie zdejmować maski bez uprzedniego odpięcia dolnych pasków uprzęży (e). 10. Σημείωση: μην αφαιρείτε τη μάσκα χωρίς να έχετε αποσυνδέσει προηγουμένως τις κάτω ζώνες του ακουστικού (e). 11. не снимайте маску, не расстегнув нижних ремней пращи(e). 12. Poznámka: masku neodstraňujte, aniž byste nejřivně odpojili dolní popruhy krytu (e). 13. 備註: 如果沒有解開頭套 (e) 的下部繫帶, 請勿取下鼻罩。 14. 注: 如需移除鼻罩, 请先松开衬垫下方的皮带(e). 15. 備考: ヘッドギア(e)の下ベルトを外さずにマスクを取り外さないでください。

**E**


4


5

6


**E**


7

8


9

# F

1. COMMENT LAVÉR LE MASQUE AU QUOTIDIEN 2. HOW TO WASH THE MASK DAILY 3. KORREKTES TÄGLICHES WASCHEN DER MASCHINE 4. CÓMO LAVAR A DIARIO LA MÁSCARA 5. COME LAVARE QUOTIDIANAMENTE LA MASCHERA 6. COMO LAVAR QUOTIDIANAMENTE A MÁSCARA 7. TVÄTTA MASKEN DAGLIGEN SÅ HÄR 8. HET MASKER DAGELIJKS WASSEN 9. JAK WYKONYWAĆ CODZIENNE CZYSZCZENIE MASKI 10. ΠΩΣ ΠΛΕΝΟΥΜΕ ΤΗ ΜΑΣΚΑ ΣΕ ΚΑΘΗΜΕΡΙΝΗ ΒΑΣΗ 11. ЕЖЕДНЕВНАЯ ОЧИСТКА МАСКИ 12. JAK MASKU DENNĚ MÝT 13. 如何日常清洗鼻罩 14. 如何日常清洗鼻罩 15. マスクの毎日の洗浄方法

**1.** Faites attention lorsque vous retirez la saleté des orifices d'exhalation du coude anti-asphyxie (c). **2.** Be very careful to remove any dirt from the exhalation holes in the anti-asphyxia fitting (c). **3.** Achten sie besonders darauf, dass die entlüftungslöcher des erstickungsschutz-anschlusses von schmutz befreit werden (c). **4.** Preste particular atención cuando elimine la suciedad de los orificios de exhalación del acoplamiento anti-asfíxia (c). **5.** Prestare particolare cura nel rimuovere lo sporco dai fori di esalazione del raccordo anti-asfissia (c). **6.** Prestar atenção especial ao remover a sujidade dos furos de exalação do encaixe anti-asfíxia (c). **7.** Var särskilt försiktig när du ska avlägsna smuts från utandningshålen på nödutandningsventilens koppling (c). **8.** Het valij zorgvuldig uit de uitademingsopeningen van de anti-verstikkingsaansluiting (c) verwijderen. **9.** Należy ostrożnie usuwać zabrudzenia z otworów wentylacyjnych złącza zapobiegającego uduszeniu(c) **10.** Δείξτε ιδιαίτερη φροντίδα σε ό,τι αφορά την αφαίρεση της βρωμιάς από τις οπές εκπομπής της αντιασφυξιογόνας ενώσεως (c) **11.** Особенно тщательно прочищайте отверстия для выхода выдыхаемого воздуха антиасфиксийного патрубка (c). **12.** Mimořádnou péči věnujte odstranění nečistot vydechovacích otvorů spoje proti udušení (c) **13.** 在去除防窒息接头 (C) 呼吸孔中的污垢时, 要格外小心。 **14.** 请仔细清洁防窒息连接件上呼气孔内的污渍。(C) **15.** 窒息防止コネクタの呼吸孔から汚れをしっかりと取り除いてください。(C)

**1. Remarques 2. Notes 3. Hinweise 4. Notas 5. Note 6. Notas 7. Anmärkningar 8. Opmerkingen 9. Uwagi 10. Σημειώσεις 11. Примечания 12. Poznámky 13. 注意 14. 注意 15. 注意事項**

▶ **1.** Le harnais (e) et le fermoir (f) doivent être lavés toutes les semaines. **2.** Wash the headgear (e) and buckle (f) every week. **3.** Das kopfteil (e) und die schnalle (f) müssen jede woche gewaschen werden. **4.** La cofia (e) y la hebilla (f) deben lavarse una vez a la semana. **5.** La cuffia (e) e la fibbia (f) devono essere lavate ogni settimana. **6.** O protetor de cabeça (e) e a fivela (f) devem ser lavados toda semana. **7.** Hättan (e) och spännet (f) måste tvättas varje vecka. **8.** De kap(e) e de gesp (f) moeten wekelijks worden gewassen. **9.** Uprząż (e) i sprzączka (f) muszą być czyszczone raz w tygodniu. **10.** Το ακουστικό (e) και η πόρπη (f) θα πρέπει να πλένονται σε εβδομαδιαία βάση. **11.** Прашу (e) и пряжку (f) следует мыть еженедельно. **12.** Kryt (e) a přezku (f) je třeba čistit jednou týdně. **13.** 頭套 (e) 和帶扣 (f) 必須每週清洗。 **14.** 衬垫(e)和卡扣(f)应每周清洗一次。 **15.** ヘッドギア(e)とバックル(f)は必ず週1回洗浄してください。

▶ **1.** Laissez les composants du masque à air sécher à l'abri de la lumière directe du soleil. **2.** Dry the mask components in the air away from direct sunlight. **3.** Lassen sie die maskenteile an der luft und geschützt vor direkter sonneneinstrahlung trocknen. **4.** Enjuagar los componentes de la máscara al aire libre protegidos de la luz directa del sol. **5.** Far asciugare i componenti della maschera all'aria al riparo dalla luce solare diretta. **6.** Deixar secar os componentes da máscara ao ar ao reparo da luz solar direta. **7.** Låt maskens komponenter torka i friska luften, skyddade från direkt solljus. **8.** De onderdelen van het masker in de lucht laten drogen, beschermd tegen rechtstreeks zonlicht. **9.** Osuszyć komponenty maski pozostawiając na powietrzu, z dala od bezpośredniego działania promieni słonecznych. **10.** Στεγνώστε τα εξαρτήματα της μάσκας στον αέρα, μακριά από απευθείας ηλιακή ακτινοβολία. **11.** Сушите компоненты маски в месте, защищенном от действия прямых солнечных лучей. **12.** Komponenty masky nechte vyschnout na vzduchu mimo dosah přímého slunečního záření. **13.** 將鼻罩的組件自然晾乾, 並避免陽光直射。 **14.** 请将鼻罩配件置于空气流通处风干, 避免阳光直射。 **15.** 直射日光的当たらない場所で、各マスク部品を乾かしてください。

▶ **1.** N'utilisez pas de solution contenant de l'alcool, des composés aromatiques, des humidifiants, des agents antibactériens, des agents hydratants ou des huiles aromatiques pour nettoyer le masque. **2.** Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life. **3.** Benutzen sie zur reinigung der

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## 1 - INTENDED USE

Respireo SOFT Facial is a reusable oro-nasal mask with calibrated exhalation orifices (Vented) designed to be used by a single patient at home or by multiple patients in a hospital setting.

Respireo SOFT Facial is intended to be used by adult patients (weighing over 30 kg), who have been prescribed non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel (positive airway pressure therapy - BiPAP).

 2 - WARNINGS

- 2.1 The mask must only be used with positive pressure generators recommended by a physician or respiratory therapist.
- 2.2. The anti-asphyxia valve fitting ensures patient safety even in case of failure of the CPAP or Bi-Level system. Do not use the mask if the valve is damaged or distorted.
- 2.3. The mask should only be used with the positive pressure equipment in operation. Exhaled air evacuation through the exhalation orifices is guaranteed if the positive pressure apparatus supplies a pressure of no less than 4 cmH<sub>2</sub>O. At a pressure lower than 4 cmH<sub>2</sub>O, the amount of air coming out from the exhalation orifices will not be enough and partial re-inhalation may occur.
- 2.4. Never cover the exhalation orifices. Before using the mask, always check that the orifices are not clogged, even partially, as this may lead to partial re-inhalation.
- 2.5. In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- 2.6. In case of additional oxygen administration, smoking or using open flames is prohibited.
- 2.7. When using oxygen and the positive pressure equipment is not in operation, turn off the oxygen dispenser to prevent accumulation resulting in a risk of fire.

- 2.8. Do not use the mask if there are signs of vomiting or nausea.
- 2.9. Keep away from light.
- 2.10. The mask should be cleaned and/or disinfected (see illustrations in section F).
- 2.11. Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- 2.12. In case of visible deterioration (cracks, tears, etc.), it is necessary to eliminate and replace the damaged mask component.
- 2.13. Do not leave the components of the mask unattended, some of these may be swallowed by children.
- 2.14. It is recommended to wash the mask before the first use.
- 2.15. Do not disassemble the anti-asphyxia fitting (c) from the mask body (b) under any circumstances, even when cleaning the mask.

**3 - MASK COMPONENTS (Section A)**

- a) cushion; b) mask body; c) anti-asphyxia fitting; d) swivel fitting; e) headgear; f) buckle.

**4 - MASK USE**

To use the mask correctly, see the illustrations in Sections B - C - D - E - G.

**5 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS**

**Notes:**

- The mask should always be regenerated before being used on a new patient.
- The headgear (e) must be replaced.

The mask can be regenerated for up to 20 cycles by one of the following processes:

- thermal disinfection at a high level;
- chemical disinfection at a high level;

Detailed mask regeneration instructions can be found on the Internet at [www.device.airliquidehealthcare.com](http://www.device.airliquidehealthcare.com) or from your retailer.

**6 - TROUBLESHOOTING**

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
• No air comes out of the exhalation orifices.	• Flow generator turned off or not working.	• Turn on the flow generator; • Adjust the flow generator.
• Difficulty breathing.	• Patient circuit not inserted correctly.	• Connect the patient circuit to the flow generator and the mask.
• Noisy mask.	• Vent ports are or partially obstructed  • Mask is not properly assembled  • Mask is not correctly positioned	• Clean the orifices using a soft bristle brush; • Disassemble the mask and then reassemble it as described in the user manual. • Reposition the mask on your face and adjust the headgear as described in the user manual.
• Air leaks between the cushion and the face.	• The mask is not properly positioned on the face. • The straps of the headgear are not properly adjusted (too tight or too loose). • The size of the mask may not be correct.	• Reposition the mask and adjust the headgear as described in the user manual.  • Contact your doctor or the technician to find the correct size


## 7 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Mask body (b)	Polycarbonate
Anti-asphyxia fitting (c)	Polycarbonate + Silicone
Swivelling coupling (d)	Co-polyester
Headgear (e)	Nylon + Polyurethane
Buckle (f)	Nylon

## 8 - TECHNICAL DATA

### PRESSURE CURVE - FLOW

**Note:** Due to machining differences, the outflow from the exhalation orifices may vary ( $\pm 10\%$ ).

<b>Standards</b>	CE marking in accordance with European Directive 93/42/EEC, Class IIa and subsequent implementation modifications - UNI EN ISO 14971, UNI EN ISO 5356-1, UNI EN ISO 10993-1, UNI EN ISO 17510, UNI EN ISO 17664
<b>Treatment pressure</b>	4 cmH <sub>2</sub> O – 40 cmH <sub>2</sub> O
<b>Dead space (size L)</b>	290 ml
<b>Resistance</b>	at 50 l/min: 0.40 cmH <sub>2</sub> O at 100 l/min: 0.70 cmH <sub>2</sub> O
<b>Resistance in single fault condition, when the anti-asphyxia valve is opened in the atmosphere</b>	Inhalation at 50 l/min: 1,5 cmH <sub>2</sub> O Exhalation at 50 l/min: 1,3 cmH <sub>2</sub> O
<b>Opening pressure of the anti-asphyxia valve</b>	1.8 cmH <sub>2</sub> O
<b>Closing pressure of the anti-asphyxia valve</b>	3.0 cmH <sub>2</sub> O
<b>Weighted sound pressure A</b>	30.9 $\pm$ 3.2 dB(A)
<b>A-Weighted Sound Power (UNI EN ISO 3744)</b>	38.9 $\pm$ 3.2 dB(A)
<b>Overall dimensions maximum (without mm (P))</b>	174 mm (A) x 116 mm (L) x 115 mm (P)
<b>Weight (size L)</b>	113 g
<b>Storage temperature</b>	-25°C / +50°C
<b>Storage humidity</b>	0% / 93%

## 9 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its

components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. The manufacturer shall be deemed responsible for the safety, reliability and performance of the device if repairs are carried out at the manufacturer or at a laboratory authorized by the manufacturer.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer. For more information on the rights under the warranty, please contact your retailer.

#### COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

#### TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.