

Domiflo

Cultivating trust



www.airliquidemedicalsystems.com



Domiflo is a pin-index pressure regulator-flow meter, designed for patients with respiratory failure. Thanks to his compactness and lightness, it is ideal for homecare applications and guarantees day-to-day patient autonomy.

Functional design

- Easy to grip
- Weight-optimised
- Reduced size
- Easy cleaning

Intuitive to use

- Straightforward installation
- Smooth locking system
- Simple flow rate selection & reading

Reassuring robustness

- Shock-resistant
- Selected components
- No preventive maintenance required



Technical specifications

Available gas	Oxygen
Flowrate (L/min)	0 - 1 - 1.5 - 2 - 2.5 - 3 - 6 - 7
Max. inlet pressure	200 bar
Outlet pressure	2 bar (+ or - 0.2 bar)
Weight	480 gr
Magnetic material	29 gr
Dimensions	91 x 108 x 59 mm
Materials	ABS (plastic), brass, silicone, EPDM, zamak
Particle filtration	> 40 micrometres
Standards compliance	NF EN ISO 10524-1 / NF EN ISO 407
Lifetime	5 years

Other versions available upon request

Reference

Domiflo CM076000



Air Liquide's Healthcare business supplies medical gases & equipment to hospitals, medical services to patients at home and is active in the fight against nosocomial infections. Its 7,800 employees serve 5,000 hospitals and 500,000 patients, on five continents. Air Liquide is the world leader in gases for industry, health and the environment.

Parc de Haute Technologie - 6, rue Georges Besse - 92182 Antony Cedex - France - Tel: +33 (0)1 40 96 66 00 - Fax: +33 (0)1 40 96 67 00