

**easy-resp™**



## ORAL MASK


*WITHOUT INTRAORAL MOUTHPIECE*

**easy-resp™** is a single-patient oral mask without deliberate leakage used for NIV in hospitals and at home.

*Patented in U.S.A and in Europe*

[www.taema.com](http://www.taema.com)

[www.airliquide.com](http://www.airliquide.com)



✓ **easy-resp™ harness is patented**

- It ensures optimal stability and comfort for the patient
- Easy to put in and take off
- Strap settings only need to be done once per patient
- Only one movement is needed to take off the mask from the patient head.
- The harness remains connected to the mask

✓ **easy-resp™ is a helpful solution for the doctor**

- Increases the ventilatory efficiency of NIV
- Reduced dead space
- Excellent airtightness allows NIV+PEEP even for patients with a beard
- Efficient NIV in the event of reduced nasal permeability, whether temporary (congestion, head cold) or permanent (deviated nasal septum)
- Compatible with the Alpha 200

✓ **easy-resp™ is more convenient for the patient**

- Completely clear field of vision
- Eliminates feelings of claustrophobia
- No contact with nose: no risk of injury to bridge of the nose or leakage in the eyes
- Comfortable, flexible, small and lightweight
- During ventilation sessions, the patient can wear glasses, read, and even watch television.

## Catalogue Item

## Characteristics

<b>easy-resp™</b> mask (with medium harness)	KM209400	<input type="checkbox"/> Made of pure silicone, guaranteed latex-free <input type="checkbox"/> Class 2A CE Marking <input type="checkbox"/> Female end fitting, 22 mm diameter, adaptable to all filters <input type="checkbox"/> Standard ventilation fitting
<b>easy-resp™</b> mask (with small harness)	KM209500	
<b>easy-resp™</b> mask (with large harness)	KM209600	
- Medium harness	KM209700	
- Small harness	KM210800	
- Large harness	KM210900	
- Swivel connector	22M/15F	

### Taema

Parc de Haute Technologie  
6, rue Georges Besse, CE 80 – 92182 Antony Cedex – FRANCE  
Tél. : +33 (0)1 40 96 66 00 – Fax : +33 (0)1 40 96 67 00

The Healthcare division within the Air Liquide Group, with its 8,000 professionals worldwide, is present in medical equipment, homecare, medical gases, medical hygiene, the pharmaceutical and cosmetic excipients. Taema is part of the medical equipment business unit of the Air Liquide group's Healthcare division.