

EVA 300

Integrated heated humidifier



EVA 300 is an integrated humidifier compatible with the **Monnal T20 and T30**. It is a solution that guarantees optimum comfort for the patient.



- Humidifier engineered to **protect** the ventilator from **water** in any positions
- Can be easily **dismantled** for regular cleaning such as in the dishwasher
- Water reservoir can be **filled** without having to dismantle the entire unit via an easily accessible filling port
- Water **capacity** of 300 ml assures humidification throughout the night
- Five heat settings that can be easily adjusted for **night-time** use
- A docking station for the **integration** of the humidifier with the entire unit

Technical specifications

Dimensions: (h x L x l) 125,4 x 393 x 196 mm

Weight: 0,8Kg

Protection: IP32

Class: IIa

References

EVA 300 (110V)	KC030900
EVA 300 (220V)	KC031000
EVA 300 reservoir	KY669100
Monnal T20	KC030500
Monnal T30	KC030600
Patient Circuit	KT600300

