

Selectaflo PR Domiflo



Home reducer- flowmeters

The home reducer-flowmeters allow a gas decompression and assure a good delivery for the patient.

The ergonomically designed flow selector suits every hand. It is also designed for easy cleaning.

Optimum safety is ensured thanks to the innovations in terms of new materials, a built-in safety valve and a reduced risk of adiabatic compression with a heat shield.



Selectaflo PR

The home reducer-flowmeter

Three ranges of calibrated-aperture reducer-flowmeters: 0 to 3 litres, 0 to 6 litres and 0 to 15 litres, with threaded or yoke connection to the cylinder.

Simple to use

Designed for easy flowrate adjustment. Pressure readings are also very easy to see in upright or reclining position and optimum operation is guaranteed no matter what the cylinder position.

Precision

Its exclusive dual pressure reduction system provides tremendous precision on delivered flowrates. This precision is maintained no matter what pressure variations there might be in the cylinder and despite the usual counterpressures due to the variety of interfaces (oxygen therapy nasal cannula or humidifier).

Part Numbers

• SELECTAFLO PR 0-3L paediatric	CM 067300
• SELECTAFLO PR 0-6L	CM 062700
• SELECTAFLO PR 0-6L yoke	CM 063100
• SELECTAFLO PR 0-15L	CM 062800
• SELECTAFLO PR 0-15L yoke	CM 063200
• SELECTAFLO PR 0-6L DIN	CM 064900
• SELECTAFLO PR 0-15L DIN	CM 065500



Domiflo

Outpatient reducer-flowmeter

Domiflo is the reducer-flowmeter with a calibrated nozzle.

Discreet and ergonomic shape

Lightweight and small reducer, especially designed for outpatient oxygen therapy.

No adjustments

The flowrate set by the output nozzle makes it simple to use and means that no handling is required.

Part Numbers

• DOMIFLO 2L yoke	CM 070000
• DOMIFLO 6L yoke	CM 071300
• Dual outlet selector (12 x 125)	DG 062800
• Dual outlet selector (9 x 16)	DG 062900

Taema

Parc de Haute Technologie
6, rue Georges Besse, CE 80 - 92182 Antony Cedex - FRANCE
Tél. : + 33 (0)1 40 96 66 00 - Fax : +33 (0)1 40 96 67 00

The Healthcare Division within the Air Liquide Group, with its 6000 professionals worldwide, is present in homecare, medical gases, medical hygiene, the pharmaceutical and cosmetic excipients