


SHUT-OFF VALVES

BS valve

Quarter-turn ball valve is used to isolate vacuum and medium pressure gas systems.



8 x 10 BS valve



Valve with lead-sealable plastic box

Performance

- Excellent seal tightness with the help of 2 seats and 2 compression type grommets.
- Double protection: chromium and nickel coated.
- Hard chromium plated brass ball.
- **BS VALVES** made of PTFE that do not stick.

Features

Maximum operating pressure: 20 bar.
Connections: brasing on copper tube.
Material: body in nickel brass and hard chromium plated spherical ball.
CE 0120 class IIa.

Maintenance

Annual sealing check.
Does not require preventive maintenance.

references

| Reference | Designation | Ø (mm) |
|-----------------|-------------|---------|
| BS001000 | BS VALVE | 8 x 10 |
| BS001200 | BS VALVE | 10 x 12 |
| BS001400 | BS VALVE | 12 x 14 |
| BS001600 | BS VALVE | 14 x 16 |
| BS001800 | BS VALVE | 16 x 18 |
| BS002000 | BS VALVE | 18 x 20 |
| BS002200 | BS VALVE | 20 x 22 |
| BS002800 | BS VALVE | 26 x 28 |
| BS003200 | BS VALVE | 30 x 32 |
| BS004200 | BS VALVE | 40 x 42 |
| BS005200 | BS VALVE | 50 x 52 |

accessories

| Reference | Designation |
|-----------------|--|
| BB010600 | Lead-sealable plastic box for BS valve up to 20 x 22 |
| YA011100 | "Neutral" compact sealed |